

Polarity Sequential IHC for Adult CNS

- All tissues and solutions are at room temperature (RT), unless noted. Always protect tissue from light exposure.
- For details on dissection and fixation see FlyLight Protocol - Adult Dissection and 2% Fixation.
- For mounting and embedding instructions refer to FlyLight Protocol – DPX Mounting.
- For videos of dissection of adult brains see Adult Brain Dissection or for adult CNS see Adult CNS dissection.
- For videos of mounting for DPX embedding of adult CNS see Adult Mounting or for larval CNS see Larval Mounting.
- For video demonstrations of DPX embedding see the movie DPX Embedding.

1. **Dissect.** Dissect adult brains or CNS in cold Schneider's Insect Medium (S2).
2. **Fix.** Transfer tissue to 2 mL Protein LoBind tubes filled with 2% paraformaldehyde (PFA) in S2 at RT. Fix for 55 minutes at RT while nutating.
3. **Post-fix wash.** Remove the fix and add 1.75 mL phosphate buffered saline with 0.5% Triton X-100 (PBT) and wash for a total of 4 X 10-minutes washes while nutating. If needed, store tissue in 0.5% PBT at 4°C while nutating or lay tube flat and rotate.
4. **Block Goat Serum (GS).** Remove PBT and add 200 µL 5% GS in PBT per tube. Incubate for 1.5 hours at RT on a rotator with tubes upright.
5. **Reference primary antibodies.** Remove block and add primary antibody diluted in 5% GS in PBT for a volume of 200 µL per tube. Incubate for 4 hours at RT on a rotator with tubes upright. Then continue incubation at 4°C on a rotator with tubes upright for 2 overnights.
 - Mouse nc82 (1:30 or 33.3 µL/mL)
6. **Post-reference primary washes.** Remove the primary antibody and do a brief rinse with 1.75 mL 0.5% PBT. Allow the tissue to settle to the bottom and then remove the rinse solution and add 1.75 mL 0.5% PBT. Wash for a total of 3 X 30-minute washes while nutating (or 4 X 15 minutes).
7. **Reference secondary antibody.** Remove PBT and add the reference primary antibody diluted in 5% GS in PBT for a volume of 200 µL per tube. Incubate for 4 hours at RT on a rotator with tubes upright. Continue incubation at 4°C on a rotator with tubes upright for 2-3 overnights.
 - Cy2 Goat α-Mouse (1:600 or 1.67 µL/mL)
8. **Second Fixation.** Transfer tissue to 2 mL Protein LoBind tubes filled with 2% paraformaldehyde (PFA) in PBS at RT. Fix for 55 minutes at RT while nutating
9. **Post-fix wash.** Remove the fix and add 1.75 mL phosphate buffered saline with 0.5% Triton X-100 (PBT) and wash for a total of 4 X 10-minutes washes while nutating. If needed, store tissue in 0.5% PBT at 4°C while nutating or lay tube flat and rotate.
10. **Block #2 GS.** Remove PBT and add 200 µL of 5% GS in PBT per tube. Incubate for 1.5 hours at RT on a rotator with tubes upright

11. **Neuron primary antibodies.** Remove block and add primary antibodies diluted in 5% GS in PBT for a volume of 200 μ L per tube. Incubate for 4 hours at RT on a rotator with tubes upright. Then continue incubation at 4°C on a rotator with tubes upright for 2 overnights.
 - Rat α -FLAG Tag (1:100 or 10 μ L/mL)
 - Rabbit α -HA Tag (1:600 or 1.67 μ L/mL)
12. **Post-neuron primary washes.** Remove the neuron primary antibodies and do a brief rinse with 1.75 mL 0.5% PBT. Allow the tissue to settle to the bottom and then remove the rinse solution and add 1.75 mL 0.5% PBT. Wash for a total of 3 X 30-minute washes while nutating (or 4 X 15 minutes).
13. **Neuron secondary antibodies.** Remove the PBT and add the secondary antibodies diluted in 5% GS in PBT for a volume of 200 μ L per tube. Incubate for 4 hours at RT on a rotator with tubes upright. Then continue incubation at 4°C on a rotator with tubes upright for 2-3 overnights.
 - ATTO647N Goat α -Rat (1:150 or 6.6 μ L/mL)
 - Cy3 Goat α -Rabbit (1:1000 or 1 μ L/mL)
14. **Post-secondary washes.** Remove the secondary antibodies and do a brief rinse with 1.75 mL 0.5% PBT. Allow the tissue to settle to the bottom and then remove the rinse solution and add 1.75 mL 0.5% PBT. Wash for a total of 3 X 30-minute washes while nutating (or 4 x 15 minutes). If needed, store tissue in 0.5% PBT at 4°C while nutating or lay tube flat and rotate.
15. **Pre-embedding fixation.** Remove PBT and add 1.75 mL 4% PFA in PBS at RT. Fix for 4 hours at RT while nutating.
16. **Post-4% PFA washes.** Remove the 4% PFA and do a brief rinse with 1.75 mL 0.5% PBT. Allow the tissue to settle to the bottom and then remove the rinse solution and add 1.75 mL 0.5% PBT. Wash for a total of 4 X 15-minute washes while nutating. If needed, store tissue in 0.5% PBT at 4°C while nutating or lay tube flat and rotate.
17. **Mount.** Mount the tissue on a poly-L-lysine (PLL) coated cover glass.
 - For making PLL see FlyLight Recipe – Poly-L-Lysine.
18. **Dehydrate.** Move the cover glass through a series of 7 cover glass staining jars filled with increasing concentrations of ethanol (30%, 50%, 75%, 95%, 100%, 100%, 100%). Soak the cover glass for 10 minutes in each jar.
19. **Xylene clearing.** (IN THE HOOD). Move the cover glass through a series of 3 jars filled with xylene. Soak the cover glass for 5 minutes in each jar.
20. **DPX embedding.** Add 7 drops of dibutyl phthalate in xylene (DPX) on top of the tissue mounted on the cover glass. Place the cover glass (DPX down) on a prepared slide with spacers. Use the edge of a glass slide to gently press down on the center of the cover glass to seat the cover glass onto the slide. Let the slide dry in the hood for 2 days before viewing.

Reporter Genotype

- pJFRC51-3xUAS-Syt::smGFP-HA in su(Hw)attP1; pJFRC225-5xUAS-IVS-myr::smGFP-FLAG in VK00005
- For a details on polarity constructs, please refer to Aso *et al.* 2014.
<http://elifesciences.org/content/3/e04577>

Reagents and Supplies

- Cy2 Goat α -Mouse. Jackson Immuno Research. # 115-225-166
- Cy3 Goat α -Rabbit. Jackson Immuno Research. # 111-165-144
- ATTO 647N Goat α -Rat IgG (H&L) Antibody. Rockland. # 612-156-120
- DPX Mountant for Microscopy. Electron Microscopy Sciences. # 13512, 500 mL
- Ethanol, ACS reagent, >99.5% (200 proof). Sigma Aldrich. # 459844-1L
- GS – Goat Serum. Life Technologies. # 16210-064, 100 mL
- nc82 – Mouse α -bruchpilot. Developmental Studies Hybridoma Bank. # nc82-s
- PBS - Phosphate Buffered Saline, 1X. Cellgro. # 21-040
- PFA – Paraformaldehyde. 20% PFA. Electron Microscopy Sciences. # 15713-S
- Poly-L-Lysine. Sigma Aldrich. # P1524-25MG
- Protein LoBind Microcentrifuge Tubes. Eppendorf. # 022431102
- Rat α -FLAG Tag (DYKDDDDK Epitope Tag). Novus Biologicals. # NBP1-06712
- Rabbit α -HA Tag. Cell Signal Technologies. # 3724S
- S2 – Schneider's Insect Medium. Sigma Aldrich. # S01416
- Triton X-100. Sigma Aldrich. # X100
- Xylenes. Fisher Scientific. # X5-500

Imaging Protocol - Polarity

Track 1 Ch 1	Cy2	498-543 nm	Neuropil (reference)
Track 1 Ch2	Cy5	654-735 nm	Neuron
Track 2 ChS1	Cy3	585-623 nm	Presynaptic Terminals
Dichromatic Mirror	MBS 488/561/633		
		20X	63X
	Resolution	1024 x 1024	1024 x 1024
	Pixel size	.52 x .52	.19 x .19
	Speed (pixel dwell)	7 (1.58 μ s)	9 (0.79 μ s)
	Bit	12	12
	Direction	Bidirectional \leftrightarrow	Bidirectional \leftrightarrow
	Average	1	1
	Zoom	0.8	0.7
	Pinhole (488)	38	68
	Interval	1 μ m	0.38 μ m