

### Polarity Sequential IHC for Adult CNS

- Overnight incubations can typically be extended to 2-3 days, along with elimination of room temperature incubation, but total procedure time should be limited.
- All tissues and solutions are at room temperature (RT), unless noted. Always protect tissue from light exposure.
- For details on dissection and fixation see FlyLight Protocol - Adult Dissection and 2% Fixation.
- For mounting and embedding instructions refer to FlyLight Protocol – DPX Mounting.
  
- For videos of dissection of adult brains see Adult Brain Dissection or for adult CNS see Adult CNS dissection.
- For videos of mounting for DPX embedding of adult CNS see Adult Mounting or for larval CNS see Larval Mounting.
- For video demonstrations of DPX embedding see the movie DPX Embedding.

1. **Dissect.** Dissect adult brains or CNS in cold Schneider's Insect Medium (S2).
2. **Fix.** Transfer tissue to 2 mL Protein LoBind tubes filled with 2% paraformaldehyde (PFA) in S2 at RT. Fix for 65 minutes at RT while nutating.
3. **Post-fix wash.** Remove the fix and add 1.75 mL phosphate buffered saline with 0.5% Triton X-100 (PBT) and wash for a total of 4 X 10-minutes washes while nutating. If needed, store tissue in 0.5% PBT at 4°C while nutating or lay tube flat and rotate.
4. **Block Goat Serum (GS).** Remove PBT and add 200 µL 5% GS in PBT per tube. Incubate for 1.5 hours at RT on a rotator with tubes upright.
5. **Reference primary antibodies.** Remove block and add primary antibody diluted in 5% GS and 0.03% sodium azide in PBT for a volume of 200 µL per tube. Incubate for 4 hours at RT on a rotator with tubes upright. Then continue incubation at 4°C on a rotator overnight.
  - Mouse nc82 (1:30 or 33.3 µL/mL)
6. **Post-reference primary washes.** Remove the primary antibody and do a brief rinse with 1.75 mL 0.5% PBT. Allow the tissue to settle to the bottom and then remove the rinse solution and add 1.75 mL 0.5% PBT. Wash for a total of 3 X 30-minute washes while nutating.
7. **Reference secondary antibody.** Remove PBT and add the reference primary antibody diluted in 5% GS and 0.03% sodium azide in PBT for a volume of 200 µL per tube. Incubate for 4 hours at RT on a rotator with tubes upright. Then continue incubation at 4°C on a rotator overnight.
  - AF568 Goat α-Mouse (1:400 or 2.5 µL/mL)
8. **Second Fixation.** Transfer tissue to 2 mL Protein LoBind tubes filled with 2% paraformaldehyde (PFA) in PBS at RT. Fix for 60 minutes at RT while nutating
9. **Post-fix wash.** Remove the fix and add 1.75 mL phosphate buffered saline with 0.5% Triton X-100 (PBT) and wash for a total of 4 X 10-minutes washes while nutating. If needed, store tissue in 0.5% PBT at 4°C while nutating or lay tube flat and rotate.
10. **Block #2 GS.** Remove PBT and add 200 µL of 5% GS in PBT per tube. Incubate for 1.5 hours at RT on a rotator with tubes upright

11. **Neuron primary antibodies.** Remove block and add primary antibodies diluted in 5% GS and 0.03% sodium azide in PBT for a volume of 200  $\mu$ L per tube. Incubate for 4 hours at RT on a rotator with tubes upright. Then continue incubation at 4°C on a rotator overnight.
  - Rabbit  $\alpha$ -GFP (1:1000 or 1  $\mu$ L/mL)
  - Rat  $\alpha$ -HA Tag (1:100 or 10  $\mu$ L/mL)
12. **Post-neuron primary washes.** Remove the neuron primary antibodies and do a brief rinse with 1.75 mL 0.5% PBT. Allow the tissue to settle to the bottom and then remove the rinse solution and add 1.75 mL 0.5% PBT. Wash for a total of 3 X 30-minute washes while nutating.
13. **Neuron secondary antibodies.** Remove the PBT and add the secondary antibodies diluted in 5% GS and 0.03% sodium azide in PBT for a volume of 200  $\mu$ L per tube. Incubate for 4 hours at RT on a rotator with tubes upright. Then continue incubation at 4°C on a rotator overnight
  - ATTO 647N Goat  $\alpha$ -Rat (1:600 or 1.67  $\mu$ L/mL)
  - AF488 Goat  $\alpha$ -Rabbit (1:800 or 1.25  $\mu$ L/mL)
14. **Post-secondary washes.** Remove the secondary antibodies and do a brief rinse with 1.75 mL 0.5% PBT. Allow the tissue to settle to the bottom and then remove the rinse solution and add 1.75 mL 0.5% PBT. Wash for a total of 3 X 30-minute washes while nutating. If needed, store tissue in 0.5% PBT at 4°C while nutating or lay tube flat and rotate.
15. **Pre-embedding fixation.** Remove PBT and add 1.75 mL 4% PFA in PBS at RT. Fix for 4 hours at RT while nutating.
16. **Post-4% PFA washes.** Remove the 4% PFA and do a brief rinse with 1.75 mL 0.5% PBT. Allow the tissue to settle to the bottom and then remove the rinse solution and add 1.75 mL 0.5% PBT. Wash for a total of 4 X 15-minute washes while nutating. If needed, store tissue in 0.5% PBT at 4°C while nutating or lay tube flat and rotate.
17. **Mount.** Mount the tissue on a poly-L-lysine (PLL) coated cover glass.
  - For making PLL see FlyLight Recipe – Poly-L-Lysine.
18. **Dehydrate.** Move the cover glass through a series of 7 cover glass staining jars filled with increasing concentrations of ethanol (30%, 50%, 75%, 95%, 100%, 100%, 100%). Soak the cover glass for 10 minutes in each jar.
19. **Xylene clearing.** (IN THE HOOD). Move the cover glass through a series of 3 jars filled with xylene. Soak the cover glass for 5 minutes in each jar.
20. **DPX embedding.** Add 7 drops of dibutyl phthalate in xylene (DPX) on top of the tissue mounted on the cover glass. Place the cover glass (DPX down) on a prepared slide with spacers. Use the edge of a glass slide to gently press down on the center of the cover glass to seat the cover glass onto the slide. Let the slide dry in the hood for 2 days before viewing.

## Reporter Genotype

- UAS-Syt-HA, 20XUAS-CsChrimson-mVenus trafficked in attP18
- UAS-Syt-HA: Robinson et al., 2002; <https://doi.org/10.1038/nature00915>
- 20XUAS-CsChrimson-mVenus trafficked in attP18: Klapoetke et al., 2014; <http://dx.doi.org/10.1038/nmeth.2836>

## Reagents and Supplies

- nc82 – Mouse  $\alpha$ -bruchpilot. Developmental Studies Hybridoma Bank. # nc82-s
- Rabbit Polyclonal  $\alpha$ -GFP. Life Technologies. #A11122
- Rat  $\alpha$ -HA. Sigma Aldrich. # 11867423001
- AF568 Goat  $\alpha$ -Mouse. Life Technologies. #A11031
- AF488 Goat  $\alpha$ -Rabbit. Invitrogen. #A11034
- ATTO 647N Goat  $\alpha$ -Rat IgG (H&L) Antibody. Rockland. # 612-156-120
- DPX Mountant for Microscopy. Electron Microscopy Sciences. # 13512, 500 mL
- Ethanol, ACS reagent, >99.5% (200 proof). Sigma Aldrich. # 459844-1L
- GS – Goat Serum. Life Technologies. # 16210-064, 100 mL
- PBS - Phosphate Buffered Saline, 1X. Cellgro. # 21-040
- PFA – Paraformaldehyde. 20% PFA. Electron Microscopy Sciences. # 15713-S
- Poly-L-Lysine. Sigma Aldrich. # P1524-25MG
- Protein LoBind Microcentrifuge Tubes. Eppendorf. # 022431102
- S2 – Schneider's Insect Medium. Sigma Aldrich. # S01416
- Triton X-100. Sigma Aldrich. # X100
- Xylenes. Fisher Scientific. # X5-500
- Sodium azide. Sigma Aldrich. # 71289-5G

## Imaging Protocol - Polarity

Track 1 Ch 1	488	498-543 nm	Neuron ( $\alpha$ -GFP)
Track 1 Ch2	633	654-735 nm	Axon Terminals ( $\alpha$ - syt::HA)
Track 2 ChS1	561	585-623 nm	Reference ( $\alpha$ -brp)
Dichromatic Mirror	MBS 488/561/633		
		<b>20X</b>	<b>63X</b>
	Resolution	1024 x 1024	1024 x 1024
	Pixel size	.52 x .52	.19 x .19
	Speed (pixel dwell)	9 (0.79 $\mu$ s)	9 (0.79 $\mu$ s)
	Bit	12	12
	Direction	Bidirectional $\leftrightarrow$	Bidirectional $\leftrightarrow$
	Average	1	1
	Zoom	0.8	0.7
	Pinhole (488)	38	68
	Interval	1 $\mu$ m	0.38 $\mu$ m